

Intracellular Cytokine Staining Protocol

Flow Cytometry Perm Buffer (10X) Cat. No. TNB-1213-L150

Other Materials Required:

- Flow Staining Buffer (1X) Cat. No. TNB-4222-L500
- Fixation Buffer Cat. No. TNB-8222-L100

Buffer and Solution Preparation

Prepare a 1X working solution of Permeabilization Buffer by diluting the 10X concentrate in distilled water prior to use.

Stimulation of Cells (optional)

A variety of methods can be used to stimulate cells to produce cytokines and investigators should determine which protocols will allow for appropriate activation of their particular sample to meet experimental needs. Tonbo Biosciences offers a **Cell Stimulation Cocktail (TNB-4975-UL100)** that can be used to elicit cytokine production and retention in a variety of cell types allowing for detection of intracellular proteins by flow cytometry. Refer to our '*Cell Stimulation for Cytokine Production Protocol*' for further details.

Experimental Procedure

1. Aliquot cell samples to tubes in a volume and at a cell concentration suitable for staining.
2. Stain cell surface antigen(s) with the recommended optimal concentration of fluorochrome labeled antibodies.
3. Wash cells in 1-2 mL Flow Staining Buffer (1X). Centrifuge at 300-400 x *g* for 5 minutes at room temperature and discard supernatant.

4. Approximately 100 μ L residual volume will generally remain in the tube. Vortex tube to dissociate the cell pellet.
5. Fix cells by adding 100 μ L of Fixation Buffer and vortex. (Final concentration approximately 2% Formaldehyde)
6. Incubate tubes in the dark at room temperature for 20-30 minutes.
7. Wash cells in Flow Staining Buffer (1X). Centrifuge at 300-400 x *g* for 5 minutes at room temperature and discard supernatant.
8. Resuspend the cell pellet in 1-2 mL of 1X Permeabilization Buffer.
9. Incubate tubes in the dark at room temperature for 5 minutes.
10. Centrifuge samples at 300-400 x *g* for 5 minutes at room temperature room temperature for 5 minutes and discard the supernatant.
11. Resuspend the cells in 100 μ L of 1X Permeabilization Buffer. Add the recommended amount of fluorochrome-labeled antibody for detection of intracellular antigen(s) to cells and incubate in the dark at room temperature for 20-60 minutes.

Note: Antibodies for intracellular staining should always be prepared in 1X Permeabilization Buffer.

12. Add 2 mL of 1X Permeabilization Buffer to each tube.
13. Centrifuge samples at 300-400 x *g* for 5 minutes at room temperature and discard the supernatant.
14. Add 2 mL of Flow Staining Buffer (1X) to each tube.
15. Centrifuge samples as in Step 13 and discard the supernatant.
16. Resuspend stained cells in an appropriate volume of Flow Staining Buffer (1X) and acquire samples on a flow cytometer.