IntraPrep Permeabilization Reagent

REF A07803 150 tests; 6 x 5 mL
2 x 0.1 mL / test

TABLE OF CONTENTS

English	2
Français (FR)	8
Deutsch (DE)	15
Italiano (IT)	22
Español (ES)	29
Português Portugal (PT-PT)	36
Dansk (DA)	43
Svenska (SV)	50
Norsk (NO)	57
Suomi (FI)	64
Ελληνικά (EL)	71
日本語 (JP)	78
中文 (ZH-CN)	85
Lietuviškai (LT)	91
Magyar (HU)	98
Polski (PL)	105
Čeština (CZ)	112
Slovenčina (SK)	119
한국어 (KO)	126
Türkçe (TR)	132
Русский (RU)	139
eesti keel (ET)	146
Hrvatski (HR)	153
Български (BG)	160
中文 (ZH-TW)	167
Română (RO)	173
Slovenščina	180
Србија (SR)	187
Latviešu (LV)	194
Українська (UK)	201
Português Brasil (PT-BR)	208
Nederlands (NL)	215
Tiếng Việt (VI)	222
Қазақша	229
APPENDIX	236
REFERENCES	237

	Reagent 1	Reagent 2
	Fixation agent	Permeability agent
Formulation	Liquid	Liquid
Active substance	Formaldehyde	Saponine
Volume	5 mL	5 mL
Number of vials	3 vials	3 vials
Volume per test	100 µL	100 µL

IntraPrep Permeabilization Reagent

REF A07803 150 tests; 6 x 5 mL
2 x 0.1 mL / test

For *In Vitro* Diagnostic Use

INTENDED USE

IntraPrep consists of two ready-to-use reagents, which induce permeability in the cytoplasmic membrane of leucocytes for the demonstration of intracellular antigenic determinants by means of monoclonal fluorescent antibodies. IntraPrep is used to prepare biological samples for analysis by flow cytometry. It has been optimized in order to minimize the non-specific staining in this type of analysis (1,2,3,4).

PRINCIPLE

As a first step, cells are fixed with reagent 1. After washing, permeability is induced with reagent 2 and remaining erythrocytes are lysed. During this stage, the cells are brought into contact with conjugated monoclonal antibodies specific for intracellular antigenic determinants. The leucocytes are then analyzed by flow cytometry.

The demonstration of surface antigenic determinants nevertheless remains possible. In this case, the specific conjugated monoclonal antibodies are incubated before fixation.

The flow cytometer measures light diffusion and the fluorescence of cells. It makes possible the localization of the population of interest within the electronic window defined on a histogram, which correlates the orthogonal diffusion of light (Side Scatter or SS) and the diffusion of narrow-angle light (Forward Scatter or FS). Other histograms combining two of the different parameters available on the cytometer can be used as supports in the gating stage depending on the application chosen by the user.

The fluorescence of the delimited cells is analyzed in order to distinguish the positively-stained events from the unstained ones. The results are expressed as a percentage of positive events in relation to all the events acquired by the gating.

INTENDED USER

This product is intended for laboratory professional use.

SAMPLES

Venous blood or bone marrow samples must be taken using sterile tubes containing an EDTA salt or ACD or heparin as the anticoagulant.

The samples should be kept at room temperature (18 – 25°C) and not shaken. The sample should be homogenized by gentle agitation prior to taking the test sample.

The samples must be analyzed within 24 hours of venipuncture.

WARNING AND PRECAUTIONS

1. Do not use the reagent beyond the expiry date.
2. Do not freeze.
3. Minimize exposure to light.
4. Avoid microbial contamination of the reagents, or false results may occur.
5. Reagent 1 contains formaldehyde. Formaldehyde is toxic and allergenic. It is thought to be a carcinogenic agent.
6. Reagent 2 contains Sodium Azide (NaN₃). It must be handled with care. Do not take internally and avoid all contact with the skin, mucosa and eyes. Moreover, in an acid medium, sodium azide can form the potentially dangerous hydrazoic acid. If it needs to be disposed of, it is recommended that the reagent be diluted in a large volume of water before pouring it into the drainage system so as to avoid the accumulation of sodium azide in metal pipes and to prevent the risk of explosion.
7. All blood samples must be considered as potentially infectious and must be handled with care (in particular: the wearing of protective gloves, gowns and goggles).
8. Never pipet by mouth, and avoid all contact of the samples with the skin, mucosa and eyes.

9. Blood tubes and disposable material used for handling should be disposed of in ad hoc containers intended for incineration.

10. Reagents and waste should be eliminated according to local requirements.

GHS HAZARD CLASSIFICATION

Reagent 1: Fixation

DANGER

H302

Harmful if swallowed.

H313

May be harmful in contact with skin

H314

Causes severe skin burns and eye damage.

H317

May cause an allergic skin reaction.

H341

Suspected of causing genetic defects.

H350

May cause cancer.

H370

Causes damage to organs.

P201

Obtain special instructions before use.

P280

Wear protective gloves, protective clothing and eye/face protection.

P303+P361+P353

IF ON SKIN (or hair): Rinse skin with water.

P305+P351+P338

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310

Immediately call a POISON CENTER or doctor/physician.

Methanol 1 - 2%

Formaldehyde 5 - 10%

Safety Data Sheet is available at beckman.com/techdocs

STORAGE AND STABILITY

IntraPrep is stored at 18 – 25°C.

Stability of closed vial: the reagent is stable for 517 days.

Stability of opened vial: the reagent is stable for 90 days.

See lot specific Certificate of Analysis at www.beckman.com.

EVIDENCE OF DETERIORATION

Any change in the physical appearance of the reagents may indicate deterioration and the reagent should not be used.

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-742-2345 (USA or Canada) or contact your local Beckman Coulter Representative.

CONTENTS

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

MATERIALS NEEDED BUT NOT SUPPLIED WITH KIT:

- Sampling tubes and material necessary for sampling.
- Automatic pipettes with disposable tips for 10, 20, 50, 100 and 500 µL.
- Plastic haemolysis tubes.

- Specific monoclonal antibodies (mAb).
- Negative controls.
- Leucocyte fixation reagent. For example : IOTest 3 Fixative Solution (Ref. A07800).
- Buffer (PBS: 0.01 M sodium phosphate; 0.145 M sodium chloride; pH 7.2).
- Centrifuge.
- Automatic agitator (Vortex type).
- Flow cytometer.

PROCEDURE

A - Intracytoplasmic staining

The number of red blood cells present in the sample should be less than 6×10^6 per μL ($6 \times 10^{12}/\text{L}$). Dilute if necessary in PBS.

The number of leucocytes present in the sample should be less than $5 \times 10^3/\mu\text{L}$ ($5 \times 10^9/\text{L}$). Dilute if necessary in PBS.

For each sample analyzed, in addition to the test tube, one control tube can be added in which the cells are mixed in the presence of the negative control, corresponding to the specific staining chosen.

1. Add 50 μL of the test sample to each tube.
2. Add 100 μL of reagent 1 to the each tube. Vortex vigorously immediately after each addition.
3. Incubate for 15 minutes at room temperature (18 – 25°C).
4. Add 4 mL of PBS to each tube.
5. Centrifuge for 5 minutes at 300 x g at room temperature.
6. Remove the supernatant by aspiration.
7. Add 100 μL of the reagent 2 to each tube. DO NOT VORTEX, leave reagent 2 to diffuse naturally into the cell pellet.
8. Incubate for 5 minutes at room temperature (18 – 25°C) WITHOUT SHAKING.
9. Shake slowly by hand for 2 to 3 seconds.
10. Add the necessary amount of mAb (specific for an intracytoplasmic antigenic determinant) to each test tube and if necessary, the appropriate amount of the negative control to each control tube.
11. Gently vortex tube by tube.
12. Incubate for 15 minutes at room temperature (18 – 25°C) protected from light.
13. Add 4 mL of PBS to each tube.
14. Centrifuge for 5 minutes at 300 x g at room temperature.
15. Remove the supernatant by aspiration and resuspend the cell pellet in 0.5 or 1 mL of IOTest 3 Fixative Solution (Ref. A07800) at its working concentration (1X). Thus fixed, the preparations can be stored between 2 and 8°C away from light for 24 hours.

B - Intracytoplasmic and membranous staining

The number of red blood cells present in the sample should be less than 6×10^6 per μL ($6 \times 10^{12}/\text{L}$). Dilute if necessary with PBS.

The number of leucocytes in the sample should be less than $5 \times 10^3/\mu\text{L}$ ($5 \times 10^9/\text{L}$). Dilute if necessary with PBS.

For each sample analyzed, in addition to the test tube, one control tube can be added in which the cells are mixed in the presence of the negative control, corresponding to the specific staining chosen.

1. Add 50 μL of the test sample to each tube.
2. Add the necessary amount of mAb (specific for a membranous antigenic determinant) to each test tube and if necessary, add the appropriate amount of negative control to each control tube.
3. Gently vortex tube by tube.
4. Incubate for 15 minutes at room temperature (18 – 25°C) protected from light.
5. Add 100 μL of reagent 1 to the each tube. Vortex vigorously immediately after each addition.
6. Incubate for 15 minutes at room temperature (18 – 25°C).
7. Add 4 mL of PBS to each tube.
8. Centrifuge for 5 minutes at 300 x g at room temperature.
9. Remove the supernatant by aspiration.
10. Add 100 μL of the reagent 2 to each tube. DO NOT VORTEX, leave reagent 2 to diffuse naturally into the cell pellet.

11. Incubate for 5 minutes at room temperature (18 – 25°C) WITHOUT SHAKING.
12. Shake slowly by hand for 2 to 3 seconds.
13. Add the necessary amount of mAb (specific for an intracytoplasmic antigenic determinant) to each test tube and if necessary, the appropriate amount of the negative control to each control tube.
14. To each control tube, add 20 µL of isotopic control.
15. Gently vortex tube by tube.
16. Incubate for 15 minutes at room temperature (18 – 25°C) protected from light.
17. Add 4 mL of PBS to each tube.
18. Centrifuge for 5 minutes at 300 x g at room temperature.
19. Remove the supernatant by aspiration and resuspend the cell pellet in 0.5 or 1 mL of IOtest 3 Fixative Solution (Ref. A07800) at its working concentration (1X). Thus fixed, the preparations can be stored at between 2 and 8°C and away from light for 24 hours.

PERFORMANCE

Performance data are obtained using the procedure described above on less than 24 hour-old blood samples previously collected in sterile tubes with EDTA salt as anticoagulant. Analysis is performed within 2 hours following immunostaining.

PRECISION

The percent positive values were determined using Whole Blood and Bone Marrow. Each sample was run 4 times, twice a day for 1 day on 2 instruments using 2 lots of IntraPrep Permeabilization Reagent. Measurements (% positive) were made on Navios flow cytometer. Analysis was conducted based on the CLSI method EP5-A2: Evaluation of Precision Performance of Quantitative Measurement Methods.

Our acceptance criteria is dependent on the number of positive events measured for each population:

- If positive event <1500, CV <15%
- If positive event >1500, CV <10%

Whole Blood EDTA:

CD45/CD14 membranous staining = Lymphocytic Purity							
Number of positive events (Mean) = 8877							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	2.29	1.82	6.47	2.25	1.65	5.82	7.09
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

CD45/CD14 membranous staining = Lymphocytic Recovery							
Number of positive events (Mean) = 8877							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	0.40	1.00	0.81	0.39	0.25	0.65	1.34
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Anti-Myeloperoxidase-FITC intracytoplasmic staining on Monocytes							
Number of positive events (Mean) = 1070							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	1.24	2.37	4.37	1.18	1.19	4.01	5.10
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Anti-Myeloperoxidase-FITC intracytoplasmic staining on Granulocytes							
Number of positive events (Mean) = 13246							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	0.09	0.03	0.16	0.04	0.03	0.12	0.16
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Whole Blood HEPARIN:

CD45/CD14 membranous staining = Lymphocytic Purity							
Number of positive events (Mean) = 6824							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	1.33	1.17	2.74	1.65	1.66	1.73	3.40
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

CD45/CD14 membranous staining = Lymphocytic Recovery							
Number of positive events (Mean) = 6824							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	2.92	2.40	3.01	0.55	0.64	0.34	3.89
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Anti-Myeloperoxidase-FITC intracytoplasmic staining on Monocytes							
Number of positive events (Mean) = 961							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	2.02	0.96	5.17	2.32	3.74	2.94	5.74
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Anti-Myeloperoxidase-FITC intracytoplasmic staining on Granulocytes							
Number of positive events (Mean) = 20742							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	1.09	1.10	2.46	1.03	1.08	1.92	2.89
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Whole Blood ACD:

CD45/CD14 membranous staining = Lymphocytic Purity							
Number of positive events (Mean) = 12355							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	1.91	1.60	3.54	1.55	1.35	2.65	4.18
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

CD45/CD14 membranous staining = Lymphocytic Recovery							
Number of positive events (Mean) = 12355							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	0.33	0.35	1.06	0.31	0.33	0.95	1.16
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Anti-Myeloperoxidase-FITC intracytoplasmic staining on Monocytes							
Number of positive events (Mean) = 1483							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	1.52	1.27	4.80	1.18	1.51	4.30	5.11
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Anti-Myeloperoxidase-FITC intracytoplasmic staining on Granulocytes							
Number of positive events (Mean) = 32715							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	0.97	0.85	2.83	0.88	0.79	2.54	3.08
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Bone Marrow EDTA:

CD45/CD14 membranous staining = Lymphocytic Purity							
Number of positive events (Mean) = 8808							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	1.16	0.89	2.98	0.72	2.32	1.46	3.19
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

CD45/CD14 membranous staining = Lymphocytic Recovery							
Number of positive events (Mean) = 8808							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	0.37	0.61	1.34	0.31	0.90	0.92	1.50
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Anti-Myeloperoxidase-FITC intracytoplasmic staining on Monocytes							
Number of positive events (Mean) = 2212							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	1.97	1.00	4.11	1.00	1.06	3.44	4.34
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Anti-Myeloperoxidase-FITC intracytoplasmic staining on Granulocytes							
Number of positive events (Mean) = 33181							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	0.11	0.29	0.37	0.10	0.13	0.33	0.48
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

LYMPHOCYTIC PURITY AND RECOVERY

Lymphocytic purity and recovery have been evaluated according to the recommendations of the CDC (5). The blood of 10 healthy donors sampled in EDTA, Heparin and ACD and 5 bone marrows were labeled with a mixture of monoclonal antibodies CD45-FITC and CD14-PE. The mean values of the recovery and purity as well as the range are given in the following tables:

Whole Blood EDTA		
Parameter	Recovery	Purity
Average	93.6	89.7
Min/Max	86.1 / 97.6	83.5 / 95.5

Whole Blood HEPARIN		
Parameter	Recovery	Purity
Average	94.3	87.6
Min/Max	84.9 / 97.7	81.5 / 94.3

Whole Blood ACD		
Parameter	Recovery	Purity
Average	94.2	94.6
Min/Max	91.1 / 98.5	87.6 / 97.7

Bone Marrow EDTA		
Parameter	Recovery	Purity
Average	80.7	71.8
Min/Max	78.3 / 85.8	60.1 / 83.4

LIMITATIONS

1. Flow cytometry may produce false results if the cytometer has not been aligned perfectly, if fluorescence leaks have not been correctly compensated for and if the regions have not been carefully positioned.
2. Certain antigenic determinants can be sensitive to formaldehyde or to saponine. Each laboratory must validate the conditions for use of monoclonal antibodies.
3. Accurate and reproducible results will be obtained as long as the procedures used are in accordance with the technical insert leaflet and compatible with good laboratory practices.
4. This reagent has been optimized so as to offer the best specific signal/non-specific signal ratio. Therefore, it is important to adhere to the ratio of reagent volume/number of leucocytes and erythrocytes in every test.
5. In the case of hypercellularity, dilute the specimen in PBS so as to obtain less than 5×10^9 leucocytes/L and less than 6×10^{12} red blood cells/L (6).
6. In certain disease states, such as severe renal failure or haemoglobinopathies, lysis of red cells may be slow, incomplete or even impossible. In this case, it is recommended to isolate mononucleated cells using a density gradient (Ficoll, for example) prior to staining (7).

See the Appendix for examples and references.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

ADDITIONAL INFORMATION

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on In vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

REVISION HISTORY

REVISION AF:	Release date : January 2021
REVISION AW:	
Updates to comply with the Beckman Coulter Global labelling Policy and according to IVD-R (EU)2017/746 requirements:	
Added sections	Intended User, Precision, Lymphocytic Recovery and Purity, Additional Information, Revision History
Added information	See sections Warning and Precautions, Storage and Stability, Limitations
Updated sections	Warning and Precautions, GHS Hazard Classification, Storage and Stability, Evidence of Deterioration, Procedure, Performance
Removed sections	Intra-Laboratory Reproducibility

Symbols Key

Glossary of Symbols is available at beckman.com/techdocs (document number B60062)

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
APPENDIX

- EXAMPLES

The graphs below are ungated biparametric representations (Side Scatter vs Fluorescence Intensity) of normal whole blood sample. Intracellular staining is with Anti-Myeloperoxidase-FITC (MPO-FITC) conjugated monoclonal antibody and with CD79a-PE conjugated monoclonal antibody.

Acquisition is performed with a Beckman Coulter Navios flow cytometer equipped with the Navios analysis software.

REFERENCES

1. Campana, D., Thompson, J.S., Amlot, P., Brown, S., Janossy, G., "The cytoplasmic expression of CD3 antigens in normal and malignant cells of the T lymphoid lineage", 1987, *J. Immunol.*, 138, 648-655.
2. Koeffler, H.P., Ranyard, J., Pertcheck, M., "Myeloperoxidase: its structure and expression during myeloid differentiation", 1985, *Blood*, 65, 484-491.
3. Picker, L.J., Singh, M.K., Zdraveski, Z., Treer, J.R., Waldrop, S.L., Bergstresser, P.R., Maino, V.C., "Direct demonstration of cytokine synthesis heterogeneity among human memory/effector T cells by flow cytometry", 1995, *Blood*, 86, 1408-1419.
4. Lees, O., "Immunophénotypage des cellules sanguines : du prélèvement à la lecture au cytomètre : les différentes méthodes", 1996, *Revue Française des Laboratoires*, 287, 33-39.
5. Nicholson, J.K.A., Hearn, T.L., Cross, G.D., White, M.D., "1997 revised guidelines for performing CD4+ T-cell determinations in persons infected with human immunodeficiency virus (HIV)", 1997, *MMWR* 46, RR-21-29.
6. H42-A2 Vol.27 No. 16. P30. Enumeration of Immunologically Defined Cell Populations by Flow Cytometry; Approved Guideline - Second Edition.
7. Ibrahim FF, Ghannam MM, Ali FM., "Effect of dialysis on erythrocyte membrane of chronically hemodialyzed patients"., 2002. *Ren Fail.*, Nov;24(6):779-90.

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Тел. +7 (495) 228 67 00, e-mail: beckman.ru@beckman.com

IMMUNOTECH SAS A Beckman Coulter Company, 130, Avenue de Lattre de Tassigny,

BP 177, 13276 Marseille Cedex 9, France, +(33) 4 91 17 27 27

www.beckman.com