

## Configuration LSR II (1/2)

Laser	Filtres		Fluorochromes	λ Emission	Relative Brightness
	BP	LP			
Bleu 488nm	530/30	505	FVS 520	521	●●●○○
			VioBright FITC	522	●●●●●
			VioBright 515	514	●●●●●
			Vio515	514	●●●○○
			BB515	513	●●●●●
			MitoTracker Green	512	●●●○○
			Alexa 488	520	●●●○○
			FITC	519	●●○○○
			Cell trace CFSE	516	●●●○○
			GFP	509	●●○○○
			YFP	527	●●●○○
	Fluo-3	528	●●●○○		
	JC1 (mono)	525	●●●○○		
575/26	550	PE	575	●●●●●	
		IP	620	●●●●○	
		DsRed	586	●●●○○	
		JC1 (poly)	595	●●●○○	
		DHE	592	●●●○○	
695/40	685	PerCP-Cy5.5	676	●●●○○	
		BB700	695	●●●○○	
		PerCP-Vio700	700	●●●●○	
		PerCP	675	●○○○○	
		7AAD	647	●●●○○	
		PE-Cy5.5	690	●●●○○	
PE-Cy5	670	●●●●●			
780/60	735	PE-Cy7	778	●●●○○	
		PE-Vio770	770	●●●○○	
Violet 405nm	440/40	V450	450	●●○○○	
		VioBlue	452	●●●○○	
FVS 450		448	●●●○○		
LIVE/DEAD Fix Violet		456	●●●●○		
Pacific Blue		456	●●○○○		
Cell trace Violet		454	●●●●○		
Alexa Fluor 405		422	●●○○○		
BV421		422	●●●●●		
DAPI		455	●●●●○		
Hoechst 33342		455	●●●●○		
525/50	505	BV510	510	●●●○○	
		VioGreen	520	●●●●○	
		FVS 510	509	●●●○○	
		Pacific Orange	551	●○○○○	
		V500	499	●●○○○	
		BV480	477	●●●●○	
		AmCyan	489	●●○○○	
		eCFP	475	●○○○○	
585/42	555	BV605	602	●●●●○	
		Super Bright 600	600	●●●●●	
		FVS 575V	572	●●●●○	

# Configuration LSR II (2/2)

