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How Cell Culture Medium Can Decrease Cell Viability During A Flow Cytometry Cell Sorting Experiment

1. Cell culture medium is not a good choice.

Most cell culture media are formulated with a CO₂-carbonate buffering system that is optimized for tissue culture carbon dioxide partial pressure, which is higher than the typical ambient carbon dioxide partial pressure. When these kinds of buffers are left in typical atmospheric conditions, **carbon dioxide will evaporate from the medium and cause the pH of the medium to rise into the alkaline range** (in phenol red-containing formulations, the color will become more purple), which expectedly can significantly affect viability. While tissue culture medium is often used as a suspension buffer without deleterious effect, **longer cell sorts and the sensitivity of the cell type can certainly exacerbate the effect of pH on viability**. Additionally, although it's important to formulate the buffer with some kind of protein, tissue culture media may contain serum at concentrations of upwards of 10-15%, which can be deposited over time on the sample line and cause flow rate irregularities during the sort. If tissue culture medium must be used, **a HEPES-based buffering system**, which is not dependent on carbon dioxide, is a preferable and safer choice.

2. Consider adding protein to your buffer.

PBS or HBSS can work well as a base for the suspension fluid, but you should add some kind of protein, typically 1-2% BSA or FBS, to keep the cells happy. Incidentally, HBSS is formulated with a CO₂-carbonate buffering system, but this system is meant to work at atmospheric/ambient CO₂ partial pressure rather than tissue-culture CO₂ partial pressure.

3. Add EDTA to prevent clumping.

EDTA in the suspension buffer, at 1-5 mM, can be a panacea for samples that tend to form aggregates and clog the instrument. This molecule interferes with cell-cell adhesion molecules by sequestering the divalent cations that these molecules require to function. Therefore, it can greatly reduce the incidence of clumping in the sample. Preventing clogs has an indirect effect on viability—clogs cause downtime during the sort, requiring cells to remain in suboptimal conditions for longer than necessary. The faster the cells are manipulated, the lower the likelihood that they will succumb to apoptosis or necrosis. If you will be using DNase, then EGTA is a better choice for chelation to prevent clumping. EDTA will chelate Mg²⁺ ions, which are necessary for DNase, while EGTA has a much lower affinity for Mg²⁺.

4. Add DNase to deal with dead cells.

Certain kinds of samples, especially those harvested from some kinds of solid tissues, will inevitably contain a high proportion of dead cells. **If this is the case, it may be advantageous to add DNase to the suspension buffer used for the sorting process.** Dead cells will release genomic DNA into the buffer, which tends to aggregate and form clumps. These clumps can trap cells, cause clogs and result in an overall drop in recovery after the sort. Recommended concentrations vary, but 200 µg/mL has been reported in Current Protocols In Cytometry to be effective. Alternatively, debris from dead cells can be removed during the preparation either through centrifugation gradients or depletion via magnetic separation.

5. Keep everything at the right temperature.

The temperature at which the sample is maintained during preparation, staining, and sorting is another important parameter to keep cells alive. The proper temperature strongly depends on the cell type and application.

While many researchers have success preparing and sorting cells at temperatures around 0°C, this may not be the ideal temperature for every cell type.

Nevertheless, there are some general benefits for keeping the cells on ice. Most importantly, a low temperature will slow down metabolic activities and will prevent biochemical or molecular changes that may affect cytometric or downstream results.

For example, some cells, under certain conditions, **may shed or internalize antibody-bound cell surface receptors, resulting in diminished staining and resolution on the cytometer.** Keeping cells cold can significantly reduce this effect.


Additionally, under the stress of preparation, staining, and the sorting process, some cells may undergo molecular changes—gene expression or otherwise—that could affect downstream experiments or analyses. **Keeping cells cold will slow down activity and possibly mitigate any transcriptional effects brought about by the experiment.**

Finally, stressed cells may begin to apoptose, lowering your recovery. By keeping your cells cold, you may slow down or mitigate this process so they can be deposited in the collection buffer before progressing further through this pathway.

When designing your flow cytometry experiment, check the literature to see if any particular suspension buffer and temperature conditions have been successful historically in other protocols for the particular cell type you are working. Buffer with protein, EDTA, and DNase is often a better choice than using your cell culture medium alone. In terms of temperature—cold, usually 0°C, is a good place to start, but temperature can be surprisingly cell-type dependent, so make sure you do your research.

10 Things Smart Scientists Do Before Sorting Cells

Cell sorting can be a scary proposition. A precious sample is introduced into a machine that pressurizes the cells to 70 PSI, moves them past one or more lasers, vibrates the stream at 90 kHz before decelerating the cells to atmospheric pressure before they hit an aqueous surface. Many cells survive this journey. But some do not.

Here are 10 things smart scientists do to improve their cell recovery :

1. They pre-coat their catch tubes.

A smart way to improve your cell recovery is to incubate your plastic tubes with a buffer solution containing protein. This will help reduce/eliminate the charge on the plastic.

Since the droplet containing the cell is charged, it can be attracted to the charge on the plastic. This results in the droplet hitting the side of the tube wall, and the cell dying as the small volume of liquid evaporates.

To prevent this, pre-coat the tube with protein/buffer to neutralize the plastic charge. Even better, make sure that **your tubes are not polystyrene.**

2. They know the catch buffer.

Cells are going to be traveling in a buffered saline. This is not very conducive for keeping cells alive for long periods of time. The good news is that you can improve your recovering by ensuring that the catch buffer has some – but not too much – protein in it. Typically only 10-50% protein in the catch buffer is sufficient.

3. They add HEPES.

If you're sorting into media, make sure the media is HEPES buffered. Buffers like RPMI are formulated to buffer in a CO₂ atmosphere (like the atmosphere found in your lab's incubator) and, as such, don't buffer well in our normal atmosphere.

4. They keep them cold (or warm).

Smart scientists know how their cells respond to temperature differences. Some cells do not like to be kept cold and will die quickly if sorted into 4 degree Celsius buffer.

5. They use a soybean trypsin inhibitor.

Sorting adherent cells adds a level of complexity to an experiment. The cells have to be disassociated to pass through the sorter, and this is often done with trypsin. The quickest and most common neutralization method is to add FBS to the cells.

Be careful of this – while it neutralizes the trypsin effectively, it also adds back all the components that cells need to re-adhere to each other. Try soybean trypsin inhibitor instead.

6. They filter.

Nothing is worse than a clogged nozzle when sorting. It adds time to the sort and reduces efficiency (and annoys the sort operator).

Just before sorting, make sure to pass the cells through an appropriate sized filter to remove clumps and debris. The smartest scientists go as far as looking at their cells under a microscope to ensure that there are no clumps prior to sorting.

7. They use a viability dye.

Make sure to include a viability dye in your staining panel. This will help eliminate dead cells.

Using a viability dye is always a smart decision.

8. They design proper antibody panels.

When trying to define a cell population, make sure you include both positive and negative markers in your antibody panel. The use of dump channels, negative markers and multiple positive markers will help ensure that the sorted cells are what they are supposed to be.

9. They ALWAYS count their cells.

Know the cell count at the time the cells are going onto the sorter – NOT from when you first began preparing them.

Since an optimal sort speed is typically $\frac{1}{4}$ the droplet generation frequency, over concentrating the cells will reduce purity at the back end. Bring some dilution buffer with you just in case the cells are too concentrated.

10. They are aware of their threshold settings.

The higher the threshold, the easier it is to visualize the specific cell population. But this doesn't eliminate the fact that the debris and junk are still present within the cell population you're visualizing. It means that the cytometer is ignoring it.

Here's the key: whatever the cytometer ignores will end up in the final sorted population.