


KIRAVIA Blue 520™ Antibodies

A Sparkling Advancement for Flow Cytometry



KIRAVIA Dyes™

Introducing the KIRAVIA Dyes™, a brand new fluorescent chemistry to advance applications in multicolor flow cytometry and beyond. KIRAVIA is a coined term meaning “sparkling” in Japanese. KIRAVIA Blue 520™ is the first in this new family of fluorophores. KIRAVIA Dyes™ employ a unique organic backbone that separates fluorophores to minimize quenching effects, thus allowing optimal and higher fluorophore to protein (F:P) ratios, surpassing what is possible with direct conjugations of single fluorophores like FITC.


- An upgrade over FITC in terms of brightness.
- A bright and spectrally clean fluorophore that only exhibits significant spillover into Spark Blue™ 550.
- Few, if any, fluors spill into KIRAVIA Blue 520™.
- Utilize this fluor on ubiquitously expressed markers or markers with varying levels of expression.

KIRAVIA Blue 520™ conjugates are twice as bright as FITC conjugates on average, offering a brighter alternative in panel building. KIRAVIA Blue 520™ can be used to detect antigens with a variety of abundance levels, but would ideally be saved for an antigen with ubiquitous or variable expression levels in complex multicolor panels (e.g. activation markers like CD25 and CD127).

View the most up-to-date list of KIRAVIA Blue 520™ products at: biolend.com/en-us/kiravia




Amazing Resolution




Anti-human CD4 (clone SK3), conjugated to FITC (red), Alexa Fluor® 488 (blue), BD Horizon™ Brilliant Blue 515 (green), or KIRAVIA Blue 520™ (purple) was used to stain human lysed whole blood.

Nuclear Staining and Titration




(Left) Anti-human CD3 (clone UCHT1) BV421™ and anti-human T-bet (clone 4B10) KIRAVIA Blue 520™ staining on human PBMCs. (Right) Titration curve with FITC or KIRAVIA Blue 520™ conjugated to anti-T-bet antibody.

Surface Staining



Human lysed whole blood was stained with anti-human CD3 APC (clone UCHT1) and anti-human CCR7 (clone G043H7) KIRAVIA Blue 520™.


Intracellular Staining




PMA/ionomycin-stimulated (6 hours) human PBMCs stained with anti-human CD3 (clone UCHT1) APC and anti-IFN-γ (clone 4S.B3) KIRAVIA Blue 520™.

KIRAVIA Blue 520™

A Brighter Alternative to FITC with Minimal Spillover




Emission spectra of KIRAVIA Blue 520™ as run on a (top) SONY ID7000™ Spectral Cell Analyzer or (bottom) 5-laser Cytek™ Aurora Spectral Cytometer.




Spillover impact of KIRAVIA Blue 520™ into detection channels of a 5-laser Cytek™ Aurora Spectral Cytometer.

Heat Stability



Anti-human CD3 (clone UCHT1) antibody conjugated to KIRAVIA Blue 520™ was incubated at the indicated temperatures over the course of 28 days. The antibodies were then used to stain human lysed whole blood cells. S/N ratio was determined as ratio of the MFI⁺/MFI⁻ signals.



Photostability



Human PBMCs were stained with anti-human CD3 (clone UCHT1) conjugated to either FITC or KIRAVIA Blue 520™. After staining, samples were fixed and stored in FluoroFix™ Buffer through the duration of the experiment. Samples were continually exposed to light or kept in the dark and read on a BD LSR Fortessa™ instrument.

Fixative Stability

KIRAVIA Blue 520™ shows utility with most BioLegend fixation buffers. Similar to some fluorophores, KIRAVIA Blue 520™ can be sensitive to the nuclear permeabilization and fixation process if used as a surface stain conjugate. While a decrease in the MFI⁺ signal was observed when staining with the FITC and KIRAVIA Blue 520™ conjugates, the positive population could still be resolved easily. The ability of a fluorophore conjugate to survive the fix/perm process will depend on the level of antigen expression and the sample type. For additional help, please contact our technical services group.



Human PBMCs were stained with anti-human CD3 (clone UCHT1) conjugated to either FITC or KIRAVIA Blue 520™ and then fixed as indicated. (Top) Fixed Fresh samples were fixed according to their respective protocols and read on a cytometer immediately. O/N samples were fixed by the indicated buffer and stored overnight in Cyto-Last™ Buffer before being read on the cytometer. Signal to Noise (S/N) was determined as ratio of the MFI⁺/MFI⁻ signals. (Bottom) Histogram overlay of the CD3 staining on PBMCs after treatment with True-Nuclear™ (black), Cyto-Fast™ (blue), or FluoroFix™ (red) as per their respective protocols.

A guide to the fixatives used in this experiment:

- FluoroFix™ is a gentler fixation PFA-based buffer.
- Cyto-Fast™ Fix/Perm Buffer Set is a specialized buffer for cytokine detection and offers improved staining over the Intracellular Staining Permeabilization Wash Buffer.
- True-Nuclear™ Transcription Factor Buffer Set is designed for optimal staining of intranuclear targets.
- True-Phos™ Perm Buffer is an alcohol-based buffer ideal for staining of phospho protein targets.
- Cyto-Last™ Buffer is formulated for long-term storage of samples while preserving signal over time.

KIRAVIA Dyes™ and KIRAVIA Blue 520™ are trademarks of SONY Corporation. BD Horizon™ and BD LSR Fortessa™ are trademarks of Becton, Dickinson and Company. Alexa Fluor® is a trademark of Life Technologies Corporation.