

Configuration FACS Aria (1/3)

Laser	Filtre		Fluorochromes	λ Emission
	BP	LP		
Violet 405nm	450/40		Cascade Blue	419
			DyLight 405	420
			Alexa Fluor 405	421
			VioBlue	455
			Pacific Blue	450
			FVS 450	450
			BV421	422
			V450	450
			eFluor450	450
			DAPI	461
			Hoechst 33342	461
			CFP	475
			BV480	478
	510/50	502	Amcyan	489
			Pacific Green	500
			eFluor 506	506
			V500	500
			FVS 510	510
			BV510	510
	560/40	545	Qdot565	565
Pacific Orange			550	
585/42	555	Qdot 585	585	
		FVS 575V	575	
		BV605	605	
		Super Bright 600	600	
670/30	635	Qdot655	655	
		BV650	650	
		SuperBright 645	645	
710/50	685	Qdot 705	705	
		SuperBright 702	702	
		BV711	711	
		BV785	785	

Configuration FACS Aria (2/3)

Laser	Filtre		Fluorochromes	λ Emission
	BP	LP		
Bleu 488nm	530/30	502	FITC	519
			Alexa 488	520
			GFP	510
			FVS 520	520
			YFP	527
			BB515	515
	695/40	655	BB700	700
			PerCP	678
			PerCP-Cy5.5	676
			PerCP-Vio700	700
			PerCP-eFluor 710	710

Laser	Filtre		Fluorochromes	λ Emission
	BP	LP		
Jaune-Vert 561nm	582/15		Cy3	566
			dTomato	581
			FVS 570	570
			RFP	588
			PE	575
			DesRed	585
	610/20	600	IP	620
			FVS 620	620
			PE-Texas Red	616
			mCherry	610
			PE-Alexa 610	610
			PE-eFluor 610	610
			PE-Vio615	615
			eFluor 615	615
			PE-CF 594 / Alexa 594	594
	670/14	630	7AAD	647
			PE-Cy5	670
			Alexa 668	668
			PE-Alexa 680	680
	710/50	685	PE-Cy5.5	675
			PE-Alexa 700	700
			PE-Vio700	700
	780/60	735	PE-Alexa750	750
			PE-Vio770	770
PE-Cy7			778	

Configuration FACS Aria (3/3)

Laser	Filtre		Fluorochromes	λ Emission
	BP	LP		
Rouge 635nm	660/20		APC	660
			FVS 660	660
			Cy5	666
			eFluor660	660
			Alexa 647	647
			DyLight 650	650
	730/45	690	Alexa 700	695
			Alexa 680	702
			FVS 700	700
			APC-R700	700
			Cy5.5	716
			APC-Cy5.5	716
	780/60	735	APC-Alexa750	774
			FVS 780	780
			APC-Cy7	785
			APC-H7	785
APC-Vio770			770	
APC-eFluor780			780	
APC/Fire 750			750	